

# Central Neutron Detector Operations Manual

S. Niccolai, IPN Orsay  
*cnd\_manual.tex – v1.2*

October 21, 2021

## **Abstract**

This document describes the Central Neutron Detector (CND) of CLAS12, and the procedures to operate it and to monitor its correct functioning.

# Contents

<b>1</b>	<b>Description of the detector</b>	<b>1</b>
<b>2</b>	<b>Read-out electronics</b>	<b>1</b>
<b>3</b>	<b>Initial operation</b>	<b>3</b>
<b>4</b>	<b>Slow controls</b>	<b>5</b>
<b>5</b>	<b>Monitoring</b>	<b>5</b>
<b>6</b>	<b>Troubleshooting</b>	<b>7</b>
<b>7</b>	<b>Maintenance</b>	<b>7</b>
<b>8</b>	<b>Responsible personnel</b>	<b>7</b>




Figure 1: Geometry of the scintillator barrel for the Central Neutron Detector. It consists of 3 radial layers each made of 48 trapezoidal scintillator paddles.

## 1 Description of the detector

The Central Neutron Detector (CND) is the outermost detector of the Central-Detector part of CLAS12. It contains the CTOF and the Central Tracker.

The CND is a barrel of plastic scintillator bars of trapezoidal shape, all with their long sides parallel to the beam direction (Fig. 1).

The light emitted by the scintillators of the CND is read only at the backward end of each bar, with an Hamamatsu R10533 photomultiplier placed in the low-field region of the solenoid, and connected to the bar by a  $\sim 1.5\text{-m}$ -long bent light guide; the front end of the bar is connected via a “u-turn” light guide to the neighboring paddle. This way, the light emitted at the front end of one scintillator is fed through its neighboring paddle and read by the PMT connected to its end (Fig. 2). Each PMT is encased in a cylindrical shielding made up by a 1-mm-thick layer of mu-metal and a 5-mm-thick layer of mild steel.

The detector is composed by 48 azimuthal segments and 3 layers in the radial direction, for a total of 144 scintillator bars, 144 PMTs, 72 u-turn light guides, 144 bent light guides (Fig. 2).

In order to operate the PMTs, high voltages (typically in the range of 1500 V) are provided to them by multi-channel CAEN SY527 power supplies (Fig. 3, left). The HV boards adopted for the CND are CAEN A734N (16 channels, 3 kV max voltage, 3 mA max current) (Fig. 3, right).

## 2 Read-out electronics

Figure 5 shows the scheme of the read-out electronics and connectors for the CND. The signal of each PMT is sent to an active splitter. The three splitter modules used for the CND were originally developed by IPN Orsay for the G0 experiment (Hall C, Jefferson Lab). Each module is an active 64-channels splitter with unity




Figure 2: Design of the Central Neutron Detector, inserted in the CLAS12 solenoid.



Figure 3: Left: Front panel of the HV power supply for the PMTs of the CND. Right: HV boards.




Figure 4: One of the CFD boards of the CND. The red and green switches are used, respectively, to select the channel and to change its threshold. The digital display should read 7-8 when the correct threshold is set.

gain so there is no loss of amplitude. The 64 SMA inputs are placed in the back panel. In the front panel there are 8 8-channels output connectors (DMCH) for the time signals and 4 16-channels output connectors (FASTBUS) for the charge signals. The charge signal is sent from the splitter to the flash-ADC (250 VXS, 16 channels/board, made and owned by JLab). The time signal from the splitter is sent to a constant fraction discriminator (CFD) GAN'ELEC FCC8, developed for the TAPS collaboration (Fig. 4). Each CFD module is an 8-channels CAMAC unit with LEMO 00 input connectors and 2x 8-pin output connectors in differential ECL. The threshold can be set for each channel individually using a manual switch, and no walk adjustment is required for the module. The discriminated time signal then goes to the TDC (CAEN VX1290A, 32 channels/board, 25 ps/channel resolution). In total, the read-out system includes 3 splitter modules, 19 CFD modules, 5 TDC boards, and 8 ADC boards.

### 3 Initial operation

The initial operation of the CND, after each prolonged downtime, should be carried out by experts.

To turn the detector on, all the crates (CAMAC and VXS) must be turned on.

The HV mainframe is turned on by turning the key (bottom left of the front panel, see Fig. 3) and then moving UP the enable switch (a light will turn on near the “ENABLED” writing). The 144 channels can then either be turned on manually or the HV slow controls in the Counting Room can be used for this purpose. The polarity for the PMTs of the CND is negative, and the typical average voltage is




Figure 5: Scheme of the electronics for the read-out of the signals of the CND.

HV SEC5 DC R1 SLO							
#	Description	Pw	Vmon	Imon	Status	Vset (V)	Iset (uA)
0	DC_SEC5_R1_SL1_S01-08	ON	1501.00	0.08	ON	1500.00	10.0
1	DC_SEC5_R1_SL1_S09-16	ON	1501.20	0.20	ON	1500.00	10.0
2	DC_SEC5_R1_SL1_S17-24	ON	1501.00	0.00	ON	1500.00	10.0
3	DC_SEC5_R1_SL1_S25-32	ON	1501.60	0.00	ON	1500.00	10.0
4	DC_SEC5_R1_SL1_S33-48	ON	1500.80	0.16	ON	1500.00	10.0
5	DC_SEC5_R1_SL1_S49-64	ON	1500.80	0.04	ON	1500.00	10.0
6	DC_SEC5_R1_SL1_S65-80	ON	1501.20	0.00	ON	1500.00	10.0
7	DC_SEC5_R1_SL1_S81-112	ON	1500.60	0.04	ON	1500.00	10.0
8	DC_SEC5_R1_SL1_F01-08	ON	749.00	0.08	ON	750.00	10.0
9	DC_SEC5_R1_SL1_F09-16	ON	748.80	0.00	ON	750.00	10.0
10	DC_SEC5_R1_SL1_F17-24	ON	748.80	0.08	ON	750.00	10.0
11	DC_SEC5_R1_SL1_F25-32	ON	749.00	0.00	ON	750.00	10.0
12	DC_SEC5_R1_SL1_F33-48	ON	749.00	0.00	ON	750.00	10.0
13	DC_SEC5_R1_SL1_F49-64	ON	748.80	0.04	ON	750.00	10.0
14	DC_SEC5_R1_SL1_F65-80	ON	748.80	0.00	ON	750.00	10.0
15	DC_SEC5_R1_SL1_F81-112	ON	749.00	0.08	ON	750.00	10.0
16	DC_SEC5_R1_SL1_G01-32	ON	526.20	0.00	ON	526.00	10.0
17	DC_SEC5_R1_SL1_G33-112	ON	525.80	0.00	ON	526.00	10.0
18	DC_SEC5_R1_SL2_S01-08	ON	1501.20	0.12	ON	1500.00	10.0
19	DC_SEC5_R1_SL2_S09-16	ON	1501.20	0.00	ON	1500.00	10.0
20	DC_SEC5_R1_SL2_S17-24	ON	1501.00	0.08	ON	1500.00	10.0
21	DC_SEC5_R1_SL2_S25-32	ON	1501.20	0.00	ON	1500.00	10.0
22	DC_SEC5_R1_SL2_S33-48	ON	1500.80	0.28	ON	1500.00	10.0

Figure 6: HV control GUI (Novice) for the CLAS12 Drift Chambers. The one for CND will look the same way.

around 1500 V.

The CFD boards must be turned on and the threshold must be set manually with a switch, until one reads a value around 7-8 on the digital display at the top of each board (Fig. 4). The red switch selects the channel, the green switch changes its threshold.

## 4 Slow controls

The slow controls for the CND are under development. They will be necessary only for the operation of the HV power supplies for the PMTs. Figures 6 and 7 show the examples of the HV control GUIs (Novice and Expert, respectively) for the CLAS12 Drift Chambers. The one of the CND will be similar to this, as the HV power supplies of the DCs are the same as those of the CND. The novice can only turn on or off the channels setting them to a fixed value of tension, while the expert can also modify the value of the tension and current, as well as the speed to ramp up and down the HV. The monitored voltages and currents are in the column called “Vmon” and “Imon”, respectively. For voltages of 1500 V, the monitored current should be around 0.1 mA for R10533 PMTs.

## 5 Monitoring

A monitoring GUI is available to check online various performances of the CND, such as charge distributions from fADCs, TDC spectra, timing resolutions, etc. Figures 8 and 9 show examples of quantities that can be monitored with the current version of the software. Once the detector is turned on, one can check the functioning of the individual channels: in the absence of beam, cosmic rates of about a couple of hundred Hz per channel are expected.

HV SEC5 DC											
#	Description	Pw	Vmon	Imon	Status	Vset (V)	Iset (uA)	Vmax (V)	Up (V/s)	Down (V/s)	
0	DC_SECS_R1_SL1_S01-08	1501.00	0.08	ON	1500.00	1500.00	10	10	1600	25 25 50	
1	DC_SECS_R1_SL1_S09-16	1501.20	0.20	ON	1500.00	1500.00	10	10	1700	25 25 50	
2	DC_SECS_R1_SL1_S17-24	1501.00	0.00	ON	1500.00	1500.00	10	10	1700	25 25 50	
3	DC_SECS_R1_SL1_S25-32	1501.60	0.00	ON	1500.00	1500.00	10	10	1700	25 25 50	
4	DC_SECS_R1_SL1_S33-48	1500.80	0.16	ON	1500.00	1500.00	10	10	1700	25 25 50	
5	DC_SECS_R1_SL1_S49-64	1500.80	0.04	ON	1500.00	1500.00	10	10	1700	25 25 50	
6	DC_SECS_R1_SL1_S65-80	1501.20	0.00	ON	1500.00	1500.00	10	10	1700	25 25 50	
7	DC_SECS_R1_SL1_S81-112	1501.60	0.04	ON	1500.00	1500.00	10	10	1700	25 25 50	
8	DC_SECS_R1_SL1_F01-08	749.80	0.08	ON	750.00	750.00	10	10	1100	25 25 50	
9	DC_SECS_R1_SL1_F09-16	748.80	0.00	ON	750.00	750.00	10	10	1300	25 25 50	
10	DC_SECS_R1_SL1_F17-24	748.80	0.08	ON	750.00	750.00	10	10	1500	25 25 50	
11	DC_SECS_R1_SL1_F35-32	749.80	0.00	ON	750.00	750.00	10	10	1500	25 25 50	
12	DC_SECS_R1_SL1_F33-48	749.80	0.00	ON	750.00	750.00	10	10	1500	25 25 50	
13	DC_SECS_R1_SL1_F49-64	748.80	0.00	ON	750.00	750.00	10	10	1500	25 25 50	
14	DC_SECS_R1_SL1_F65-80	748.80	0.00	ON	750.00	750.00	10	10	1500	25 25 50	
15	DC_SECS_R1_SL1_F81-112	749.80	0.04	ON	750.00	750.00	10	10	1500	25 25 50	
16	DC_SECS_R1_SL1_G01-32	526.20	0.00	ON	526.00	526.00	10	10	1600	25 25 50	
17	DC_SECS_R1_SL1_G03-112	525.80	0.00	ON	526.00	526.00	10	10	600	25 25 50	
18	DC_SECS_R1_SL2_S02-08	1501.20	0.12	ON	1500.00	1500.00	10	10	1700	25 25 50	
19	DC_SECS_R1_SL2_S09-16	1501.20	0.00	ON	1500.00	1500.00	10	10	1700	25 25 50	
20	DC_SECS_R1_SL2_S17-24	1501.00	0.08	ON	1500.00	1500.00	10	10	1700	25 25 50	
21	DC_SECS_R1_SL2_S25-32	1501.20	0.04	ON	1500.00	1500.00	10	10	1700	25 25 50	
22	DC_SECS_R1_SL2_S33-48	1500.80	0.28	ON	1500.00	1500.00	10	10	1700	25 25 50	

Figure 7: HV control GUI (Expert) for the CLAS12 Drift Chambers. The one for CND will look the same way.




Figure 8: Charge from the fADC, for the two paddles (left and right) of one layer of the CND, as obtained on cosmic-rays data using the CND monitoring GUI.




Figure 9: Plots of time differences between CND paddles used to compute the time resolution, as obtained on cosmic-rays data using the CND monitoring GUI.

## 6 Troubleshooting

The HV slow control GUI has alarms that signal either trips (red alarm) or if a channel readback is far from the set value. In case an HV alarm appears for one channel of the CND, the novice can try turning the HV off and on. The current readback should go back to the nominal value. Then he should check the rates on the monitoring GUI. If the problem persists, call the expert.

## 7 Maintenance

Given their location outside of the solenoid, all the PMTs of the CND are accessible once all the CD detectors are installed. If needed, in case of malfunctioning, the PMTs can be replaced without having to remove other CD detectors.

## 8 Responsible personnel

Individuals responsible for the Central Neutron Detector are listed in Table 8.

Name	Dept.	Phone	Email	Comments
Expert on call		757-XXX-XXXX		
S. Niccolai	Orsay	+33 6 24 81 67 78	silvia@jlab.org	First contact
D. Sokhan	Glasgow	+44 7949 175725	daria@jlab.org	Second contact
D. Carman	JLab	757-269-5586	carman@jlab.org	JLab contact
S. Boyarinov	JLab	757-269-5795	boyarinov@jlab.org	DAQ
N. Baltzell	JLab	757-269-5902	baltzell@jlab.org	Slow Controls

Table 1: Personnel responsible for the Central Neutron Detector.